



UAM Motor Drive and Drum Hoop Blender HM-5

Drum Hoop Blender HM-5

The Drum Hoop Blender (Turbola Mixer) HM-5 has been designed for mixing such as powder and granules. The Drum Hoop Mixer is equally suited to research as well as industrial use in the pharmaceutical, food, chemical and plastics fields. Drum Hoop Mixer comprises of two parts.

1. The Drive Unit with drum hoop support stands. This part is connected to the UAM Motor Drive.
2. The Acrylic mixing drum with rotating fixture.

Technical Data:

Gross Capacity	: 5 lit.
Working Capacity	: 2.5 – 3 lit.
Dimension	: 500 x 300 x 270 mm
Net Weight	: 6.5 Kg.
Gross Weight	: 33 Kg.
Packing	: 680 x 430 x 430 mm

We reserve the right to make technical changes without any prior notice