

Ready to Start a Sequence in Only 30 Minutes:



Placing the samples



Scanning the sample IDs or manual input



Insert the ready-to-use columns in the carousel and press „Start“.

DEXTech 16 processes 16 samples one after another without any manual working step in between separating the PCBs and dioxins /dl-PCBs.

Sophisticated but Simple Software

- Plausibility driven programming
- Permanent overview of process status
- Intelligent solvent management
- Prioritization of samples possible
- Non-manipulable report generation
- User authorisation management

Unbeatable Safety Features

- Password protected software
- Controlled and monitored robotic operations
- Overpressure detection
- Waste sensor
- Leakage detection

You talked to your LCTech specialist:

E-mail:

.....

.....@LCTech.de



The information contained in this brochure is based on our current knowledge and has been carefully checked. However, since we continually work on the further development of our products, please accept texts, pictures and numbers on these pages as non-binding and exemplary only.

LCTech GmbH
Daimlerstraße 4
84419 Obertaufkirchen
Germany

Tel: +49 8082 2717-0
Fax: +49 8082 2717-100
info@LCTech.de
www.LCTech.de