

Hei-CHILL 1200



With rollers for placement under a lab bench and manometer for pressure display.
Comes with a cooling capacity of 1200 W*

- Clearly arranged control elements, large LED display, membrane keypad, and window for monitoring the tempering liquid level
- Automatic start and shut down feature
- Temperature range from -10 °C to +40 °C
- Temperature stability of ± 0.5 K
- Cooling capacity at +20 °C; 1200 W
- Footprint of chiller: W 450 x H 550 x D 650 mm

* Cooling output at 20 °C measured with ethanol

Hei-CHILL 1200 - Technical Data

Dimensions (w/h/d)	450 x 550 x 650 mm
Temperature range	-10 – +40 °C
Coolant volume	7 – 14 l
Cooling capacity at +20 °C (20 °C ambient temp.)	1,200 W
Temperature accuracy	±0.5K
Pump pressure max.	1.3 bar
Pump flow max.	35 l/min