

## Hei-CHILL 600



With rollers for placement under a lab bench and manometer for pressure display.  
Comes with a cooling capacity of 600 W\*

- Clearly arranged control elements, large LED display, membrane keypad, and window for monitoring the tempering liquid level
- Automatic start and shut down feature
- Temperature range from -10 °C to +40 °C
- Temperature stability of  $\pm 0.5$  K
- Cooling capacity at +20 °C; 600 W
- Footprint of chiller: W 350 x H 480 x D 595 mm

\* Cooling output at 20 °C measured with ethanol

## Hei-CHILL 600 - Technical Data

Dimensions (w/h/d)	350 x 480 x 595 mm
Temperature range	-10 – +40 °C
Coolant volume	4 – 8 l
Cooling capacity at +20 °C (20 °C ambient temp.)	600 W
Temperature accuracy	±0.5K
Pump pressure max.	1.3 bar
Pump flow max.	35 l/min