

# Celase<sup>®</sup> GMP

## Enzyme

Distributed by:



# Celase® GMP Enzyme

## Convenience In Your Lab

- A single vial containing both collagenase and a neutral protease can digest up to 280 g of adipose tissue
- Shelf life of up to 48 months

## Clarity In Your Approach

- Protocols are available from Cytori for dissociating canine, equine, human, ovine, porcine, rabbit and rodent adipose tissue
- Technical dossier is available from Cytori to ease the transition from research to clinical applications

## Confidence In Your Result

- Used in U.S. FDA approved clinical trials for chronic heart failure, hamstring injuries, osteoarthritis of the knee, and hand manifestations of scleroderma
- Produced using avian and mammalian tissue free raw materials, aseptic processes and sterile filtration under GMP guidelines to assure the lowest levels of impurities



## Next Level Research

Celase® GMP is a proprietary, blended proteolytic enzyme designed for efficient, gentle and reproducible *in vitro* dissociation of nucleated cells from adipose tissue.

## Specifications

PARAMETER	RANGE / VALUE
Total Protein	34.4 – 51.6 mg/vial
Endotoxin	≤ 50 EU/mg
Stability	48 months at -15 to -25° C
Appearance	White lyophilizate

## Ordering Information

**Catalog Number**  
1235-01

**Pack Size**  
1 vial, 35 mg



Distributed by:

**Worthington Biochemical**

**Serlabo Technologies**

[info@serlabo.fr](mailto:info@serlabo.fr)

Tél: 04.90.23.77.20



©2015 Cytori Therapeutics, Inc. All rights reserved.  
Cytori, the Cytori logo, Celase, and the Celase logo are trademarks or registered trademarks of Cytori Therapeutics, Inc. in the United States and other select countries.

ENZ-009-LIT-GLB\_A-0415