

# Intelligent Evaporation



MADE IN  
GERMANY

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a worthwhile investment.



Ground-free condenser:  
Ventilation cap with  
PTFE inlay and threaded  
connection **prevents  
stuck glass**

**Up to 40 % higher evaporation rate** with a cooling  
surface of 2,200 cm<sup>2</sup> for all XL glass condensers

Highly resistant and durable PTFE vacuum seals  
**for many years of continuous operation**

The right operating concept  
for all requirements:  
**with digital display or  
innovative touch display.**  
Also available as Control  
model with Control-Box  
for additional control of  
vacuum and Hei-CHILL  
recirculation chiller

**Patented connection system**  
with Easy-Clip to gently  
screw-off flasks and **vapor  
tube with clamping sleeve**  
prevent stuck glass parts

**The 210 °C safety heating  
bath** with ergonomic handles  
and pour spout shuts down  
when running dry or at  
over temperature > 5 °C

For quick access, all models  
also have separate knobs  
with **LED rings showing  
the active functions**

Cable connections with  
**protection class IP 67**  
protect against corrosion  
and short-circuits

**For unlimited evaporation:** In conjunction with  
the Distimatic automatic module, you can distill  
indefinitely and automatically in 24 / 7 operation. The  
affordable alternative to large rotary evaporators.

**With standby function** and  
residual heat warning at heating  
bath temperatures above 50 °C

Metal connection to the  
device console **prevents  
tilting of the bath**