

Safe

Heating and Mixing



MADE IN
GERMANY

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a worthwhile investment.

The standard **RS 232 interface** of the MR Hei-Connect and MR Hei-End **enables the documentation and controlling via PC**. A **free software is available** from our website

The extended heating capacity of **800 W** **reduces heat-up times by 35 %** compared to other units at 600 W

The chemically resistant Kera-Disk® hot plate allows for **immediate heat transfer** resulting in quick heat-up times

Hermetically-**sealed housing** **protects** all mechanical and electronic components **from aggressive environments**

In case the heating function fails, stirring will not be discontinued to **prevent bumping**



The temperature sensor consists of **two independent safety circuits** which switches off heating in case of any overtemperature situation

In case of a short-circuit, a damage or removal of the temperature sensor from the media, the **unit powers off completely**

Damage to the stirrer is categorically ruled out even if exposed to highest temperatures – all models come with **a fire-resistant aluminum die-cast housing**

A separate on/off button for heating **prevents unintentional heat-up** – the button is illuminated for visual control. In addition, the residual heat indicator prevents you from potential burning hazards